

III


10mm/mV 25mm/sec 25Hz


Pat. ID:.....

aVR

BPL CARDIART 6108T

CARDIART






Electrocardiogram (ECG) recording on grid paper. The tracing shows two leads: aVL and aVF. The aVL lead shows a small rS pattern, and the aVF lead shows a small rS pattern. The paper includes technical specifications: 10mm/mV, 25mm/sec, and 25Hz. The BPL CARDIART 61081 logo is visible.

CARDIART

Pat. ID .....

CARDIART



V1

10mm/mv 25mm/sec 25HZ


V2

BPL CARDIART 6108T

BPL

CARDIART

Pat. ID.....



V3

BPL 10mm/mv 25mm/sec  $\approx$  25HZ

V4

BPL CARDIART 61081

BPL

V5 10mm/mV 25mm/sec

Pat. ID.....

CARDIART

BPL

V6 BPL CARDIART 6108T

CARDIART

